

Product Information

MemDX™ Antibody Discovery - Mouse EGFR (25-647) Membrane Protein, Partial, -His tag

Cat. No.: **MP0047F**

This product is for research use only and is not intended for diagnostic use.

This membrane protein is Mouse EGFR (25-647). It has been tested in SDS-PAGE, ELISA. We provide this protein to facilitate your membrane protein antibody discovery and development.

Product Specifications

Host Species

Mouse

Target Protein

EGFR

Protein Length

ECD

Molecular Weight

The protein has a calculated MW of 71.2 kDa. The protein migrates as 86-110 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Sequence

AA Leu 25 - Ser 647 (Accession # Q01279-1).

Product Description

Activity

Yes

Application

SDS-PAGE, ELISA

Expression Systems

HEK293

Tag

His tag at the C-terminus

Protein Format

Soluble

Form

LYOPH

Reconstitution

Please see Certificate of Analysis for specific instructions.

Endotoxin

<1.0 EU/μg by the LAL method

Purity

>95% as determined by SDS-PAGE.

Buffer

Lyophilized from 0.22 μm filtered solution in PBS, pH7.4. Normally trehalose is added as protectant before lyophilization.

Storage

Stored at lyophilized form at -20°C or lower. Avoid repeated freeze-thaw cycles.

The antigen can be stable for 12 months in lyophilized form after storage at -20°C to -80°C, 3 months under sterile conditions after reconstitution after storage at -80°C.

Target**Target Protein**

EGFR

Full Name

epidermal growth factor receptor

Introduction

The protein encoded by this gene is a transmembrane glycoprotein that is a member of the protein kinase superfamily. This protein is a receptor for members of the epidermal growth factor family. EGFR is a cell surface protein that binds to epidermal growth factor, thus inducing receptor dimerization and tyrosine autophosphorylation leading to cell proliferation. Mutations in this gene are associated with lung cancer. EGFR is a component of the cytokine storm which contributes to a severe form of Coronavirus Disease 2019 (COVID-19) resulting from infection with severe acute respiratory syndrome coronavirus-2 (SARS-CoV-2).

Alternative Names

Wa; Erb; Err; Wa5; wa2; Erbb; Errp; wa-2; Errb1; A1552599; 9030024J15Rik;avian erythroblastic leukemia viral (v-erb-b) oncogene homolog;waved 2

Gene ID

[13649](#)

UniProt ID

[Q01279](#)