

Product Information

MemDX™ Membrane Protein Human TMEM190 (Transmembrane protein 190)

Cat. No.: **MP0050J**

This product is for research use only and is not intended for diagnostic use.

This product is a 19.3 kDa Human TMEM190 membrane protein expressed in HEK293T. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM190

Protein Length

Full-length

Protein Class

Transmembrane

Molecular Weight

19.3 kDa

TMD

1

Sequence

MLGCGIPALGLLLLLQGSADGNGIQGFFYPWSCEGDIWDRESCGGQAAIDSPNLCLRLRCCYRNGVCYHQ
RPDENVRRKHMWALVWTCSGLLLLSCSICLFWWAKRRDVLHMPGFLAGPCDMSKSVSLLSKHRGTTKTPS
TGSVPVALSKESRDVEGGTEGEGTEEGEETEGEEED

Product Description

Expression Systems

HEK293T

Tag

C-Myc/DDK

Form

Liquid

Purification

Anti-DDK affinity column followed by conventional chromatography steps

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer

25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

TMEM190

Full Name

Transmembrane protein 190

Introduction

The TMEM190 gene is conserved in chimpanzee, Rhesus monkey, dog, cow, mouse, and rat.

Alternative Names

MDAC1

Gene ID

[147744](#)

UniProt ID

[Q8WZ59](#)