

Product Information

MemDX™ Membrane Protein Human HLA-B (Major histocompatibility complex, class I, B)

Cat. No.: **MP0139J**

This product is for research use only and is not intended for diagnostic use.

This product is a 37.8 kDa Human HLA-B membrane protein expressed in HEK293T. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

HLA-B

Protein Length

Full-length

Protein Class

Transmembrane

Molecular Weight

37.8 kDa

TMD

1

Sequence

MLVMAPRTVLLLLSAALALTETWAGSHSMRYFDTAMSRPGRGEPFISVGYYDDTQFVRFDSDAASPREE
PRAPWIEQEGPEYWDRNTQIFKTNQTDRESLRNLRGYNNQSEAGSHTLQSMYGCDVGPDGRLLRGHNQY
AYDGKDYIALNEDLRSWTAADTAAQITQRKWEAARVAEQDRAYLEGTCVEWLRRYLENGKDTLERADPPK
THVTHHPISDHEATLRCWALGFYPAEITLTWQRDGEDQTQDTELVETRPAGDRTFQKWAAVVVPSEEQR
YTCHVQHEGLPKPLTLRWEPSSQSTVPIVGIVAGLAVLAVVVIGAVVAVMCRRKSSGGKGGSYSQAACS
DSAQGSQSDVSLTA

Product Description

Expression Systems

HEK293T

Tag

C-Myc/DDK

Form

Liquid

Purification

Anti-DDK affinity column followed by conventional chromatography steps

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer

25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

HLA-B

Full Name

Major histocompatibility complex, class I, B

Introduction

HLA-B belongs to the HLA class I heavy chain paralogues. This class I molecule is a heterodimer consisting of a heavy chain and a light chain (beta-2 microglobulin). The heavy chain is anchored in the membrane. Class I molecules play a central role in the immune system by presenting peptides derived from the endoplasmic reticulum lumen. They are expressed in nearly all cells. The heavy chain is approximately 45 kDa and its gene contains 8 exons. Exon 1 encodes the leader peptide, exon 2 and 3 encode the alpha1 and alpha2 domains, which both bind the peptide, exon 4 encodes the alpha3 domain, exon 5 encodes the transmembrane region and exons 6 and 7 encode the cytoplasmic tail. Polymorphisms within exon 2 and exon 3 are responsible for the peptide binding specificity of each class one molecule. Typing for these polymorphisms is routinely done for bone marrow and kidney transplantation. Hundreds of HLA-B alleles have been described.

Alternative Names

AS; B-4901; HLAB; SPDA1

Gene ID

[3106](#)

UniProt ID

[P01889](#)