

Product Information

MemDX™ Membrane Protein Human CLDN1 (Claudin 1)

Cat. No.: **MP0209J**

This product is for research use only and is not intended for diagnostic use.

This product is a 22.6 kDa Human CLDN1 membrane protein expressed in HEK293T. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

CLDN1

Protein Length

Full-length

Protein Class

Transmembrane

Molecular Weight

22.6 kDa

TMD

4

Sequence

MANAGLQLLGFI~~LA~~FLGWIGAIVSTALPQWRIYSYAGDNIVTAQAMYEG~~LW~~MSCV~~SQ~~STGQIQCKVFDSL
LNL~~SST~~LQATRALMVVGILLGVIAIFVATVGMKCMKCLE~~DDEV~~QKMRMAVIGGAIFLLAGLAILVATAWY
GNRIVQEFYDPMTPVNARYEFGQALFTGWAAASLCLLGGALLCCSCPRKTT~~SY~~TPRPPYPKPAPSSGKDY
V

Product Description

Expression Systems

HEK293T

Tag

C-Myc/DDK

Form

Liquid

Purification

Anti-DDK affinity column followed by conventional chromatography steps

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer

25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

CLDN1

Full Name

Claudin 1

Introduction

Tight junctions represent one mode of cell-to-cell adhesion in epithelial or endothelial cell sheets, forming continuous seals around cells and serving as a physical barrier to prevent solutes and water from passing freely through the paracellular space. These junctions are comprised of sets of continuous networking strands in the outwardly facing cytoplasmic leaflet, with complementary grooves in the inwardly facing extracytoplasmic leaflet. The protein encoded by this gene, a member of the claudin family, is an integral membrane protein and a component of tight junction strands. Loss of function mutations result in neonatal ichthyosis-sclerosing cholangitis syndrome.

Alternative Names

CLD1; SEMP1; ILVASC

Gene ID

[9076](#)

UniProt ID

[O95832](#)