

Product Information

MemDX™ Membrane Protein Human CH25H (Cholesterol 25-hydroxylase)

Cat. No.: **MP0239J**

This product is for research use only and is not intended for diagnostic use.

This product is a 31.6 kDa Human CH25H membrane protein expressed in HEK293T. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

CH25H

Protein Length

Full-length

Protein Class

Transmembrane

Molecular Weight

31.6 kDa

TMD

3

Sequence

MSCHNCSDPQVLCSSGQLFLQPLWDHLRSWEALLQSPFFPVIFSITTYVGFCLPFVVDILCSWVPALRR
YKIHPDFSPSAQQLPCLGQTLYQHVMFVFPVTLHWRSPALLPHEAPELLLLHHILFCLLLFDMEFF
VWHLHHKVPWLRYTFHKVHHQNSSSFALATQYMSVWELFSLGFFDMMNVTLGCHPLTTLTFHVVNIWL
SVEDHSGYNFPWSTHRLVPFGWYGGVVHDLHSHFNCNFAPYFTHWDKILGTLRTASVPR

Product Description

Expression Systems

HEK293T

Tag

C-Myc/DDK

Form

Liquid

Purification

Anti-DDK affinity column followed by conventional chromatography steps

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer

25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

CH25H

Full Name

Cholesterol 25-hydroxylase

Introduction

This is an intronless gene that is involved in cholesterol and lipid metabolism. The encoded protein is a membrane protein and contains clusters of histidine residues essential for catalytic activity. Unlike most other sterol hydroxylases, this enzyme is a member of a small family of enzymes that utilize diiron cofactors to catalyze the hydroxylation of hydrophobic substrates.

Alternative Names

C25H; cholesterol 25-hydroxylase; h25OH; cholesterol 25-monooxygenase

Gene ID

[9023](#)

UniProt ID

[O95992](#)