

Product Information

MemDX™ Membrane Protein Human TEX28 (Testis expressed 28)

Cat. No.: **MP0250J**

This product is for research use only and is not intended for diagnostic use.

This product is a 46 kDa Human TEX28 membrane protein expressed in HEK293T. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TEX28

Protein Length

Full-length

Protein Class

Transmembrane

Molecular Weight

46 kDa

TMD

1

Sequence

MVLKAEHTRSPSATLPSNVPSCRSLSSESDGPSGPSSLADGGLAHNLQDSVRHRILYLSEQLRVEKASRD
GNTVSYLKLVSADRHQVPHIQQAFQKVNQRASATIAQIEHRLHQCHQQQLQELEGCRPEGLLLMAESDP
ANCEPPSEKALLSEPPEPGGEDGPVNLPHASRPFILESRFQSLQQGTCTEDVAQQQNLLLQKVKAEELE
EAKRFHISLQESYHSLKERSLTDLQLLLESLEEKCRQALMEEQVNGRLQGQLNEIYNLKHNLACSEERM
AYLSYERAKEIWEITETFKSRISKLEMLQQVTQLEAAEHLQSRPPQMLFKFLSPRLSLATVLLVFVSTLC
ACPSSLISSRLCTCTMLMLIGLVLAQWRRAIPATDWQEWVPSRCRLYSKDSGPPADGP

Product Description

Expression Systems

HEK293T

Tag

C-Myc/DDK

Form

Liquid

Purification

Anti-DDK affinity column followed by conventional chromatography steps

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer

25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

TEX28

Full Name

Testis expressed 28

Introduction

CThe TEX28 gene is conserved in chimpanzee, Rhesus monkey, dog, cow, mouse, and rat.

Alternative Names

tTEX; MRX99; CXorf2; TEX28P1; TEX28P2

Gene ID

[1527](#)

UniProt ID

[O15482](#)