

Product Information

MemDX™ Membrane Protein Human CD99L2 (CD99 molecule like 2)

Cat. No.: **MP0279J**

This product is for research use only and is not intended for diagnostic use.

This product is a 27.8 kDa Human CD99L2 membrane protein expressed in HEK293T. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

CD99L2

Protein Length

Full-length

Protein Class

Transmembrane

Molecular Weight

27.8 kDa

TMD

1

Sequence

MVAWRS AFLVCLAFSLATLVQRGSGDFDDFNLEDAVKETSSVKQPWDHTTTTTTNRPGTTRAPAKPPGSG
LDLADALDDQDDGRRKPGIGGRERWNHVTTTTTKRPVTTTRAPANTLGNDFDLADALDDRNDNRDDGRRKPIA
GGGGFSDKDLEDIVGGGEYKPDKGKG DGRYGSND DPGSGMVAEPGTIAGVASALAMALIGAVSSYISYQQ
KKFCFSIQQGLNADYVKGENLEAVVCEEPQVKYSTLHTQSAEPPPPPEPARI

Product Description

Expression Systems

HEK293T

Tag

C-Myc/DDK

Form

Liquid

Purification

Anti-DDK affinity column followed by conventional chromatography steps

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer

25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

CD99L2

Full Name

CD99 molecule like 2

Introduction

This gene encodes a cell-surface protein that is similar to CD99. A similar protein in mouse functions as an adhesion molecule during leukocyte extravasation. Alternate splicing results in multiple transcript variants.

Alternative Names

CD99B; MIC2L1; MIC2 like 1; MIC2-like protein 1

Gene ID

[83692](#)

UniProt ID

[Q8TCZ2](#)