

Product Information

MemDX™ Membrane Protein Human MXRA8 (Matrix remodeling associated 8)

Cat. No.: **MP0289J**

This product is for research use only and is not intended for diagnostic use.

This product is a 49 kDa Human MXRA8 membrane protein expressed in HEK293T. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

MXRA8

Protein Length

Full-length

Protein Class

Transmembrane

Molecular Weight

49 kDa

TMD

1

Sequence

MALPSRILLWKLVLQSSAVLLHSGSSVPAAAGSSVSESAVSWEAGARAVLRCQSPRMVWTQDRLHDRQ
RVLHWDLRGPGGGPARRLDLYSAGEQRVYEARDRGRLELSASAFDDGNFSLIRAVEETDAGLYTCNLH
HHYCHLYESLAVRLEVTGPPATPAYWDGEKEVLAVARGAPALLTCVNRGHVWTDHRHVEEAQQVVHWDQR
PPGVPHDRADRLDLYASGERRAYGPLFLRDRVAVGADAFERGFSLRIEPLVADEGTYSCHLHHHYCG
LHERRVFHLTVAEPHAEPPIRGSPGNGSSHSGAPGPDPTLARGHNVINVIVPESRAHFFQQLGYVLATLL
LFILLVTVLLAARRRRGGYEYSDQKSGKSKGKDVNLAFAVAAGDQMLYRSEDIQLDYKNILKERAEL
AHSPLPAKYIDLKGRKENCK

Product Description

Expression Systems

HEK293T

Tag

C-Myc/DDK

Form

Liquid

Purification

Anti-DDK affinity column followed by conventional chromatography steps

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer

25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Storage

Store at +4°C for up to one week or several months at -80°C

Target

Target Protein

MXRA8

Full Name

Matrix remodeling associated 8

Introduction

Transmembrane protein which can modulate activity of various signaling pathways, probably via binding to integrin ITGAV:ITGB3 (PubMed:22492581, PubMed:23386276). Mediates heterophilic cell-cell interactions *in vitro* (By similarity). Inhibits osteoclastogenesis downstream of TNFSF11/RANKL and CSF1, where it may function by attenuating signaling via integrin ITGB3 and MAP kinase p38 (By similarity). Plays a role in cartilage formation where it promotes proliferation and maturation of growth plate chondrocytes (By similarity). Stimulates formation of primary cilia in chondrocytes (By similarity). Enhances expression of genes involved in the hedgehog signaling pathway in chondrocytes, including the hedgehog signaling molecule IHH; may also promote signaling via the PTHLH/PTHrP pathway (By similarity). Plays a role in angiogenesis where it suppresses migration of endothelial cells and also promotes their apoptosis (PubMed:23386276). Inhibits VEGF-induced activation of AKT and p38 MAP kinase in endothelial cells (PubMed:23386276). Also inhibits VTN (vitronectin)-mediated integrin ITGAV:ITGB3 signaling and activation of PTK2/FAK (PubMed:23386276). May play a role in the maturation and maintenance of the blood-brain barrier (By similarity).

Alternative Names

ASP3; matrix-remodeling-associated protein 8; limitrin

Gene ID

[54587](#)

UniProt ID

[Q9BRK3](#)