

Product Information

MemDX™ Membrane Protein Human FCAMR (Fc fragment of IgA and IgM receptor) for

Antibody Discovery

Cat. No.: **MP0354X**

This product is for research use only and is not intended for diagnostic use.

This product is a 29.15 kDa Human FCAMR membrane protein expressed in *in vitro* wheat germ expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

FCAMR

Protein Length

Full-length

Molecular Weight

29.15 kDa

TMD

1

Sequence

MDGEATVKPGEQKEVRRGREVDYSRLIAGTLPQSHVTSRRAGWKMPFLILCLLQGSSFALPQKRPHPRWLWEGSLPSRTHLRA

Product Description

Application

Antibody Production

Expression Systems

in vitro wheat germ expression system

Tag

NO

Protein Format

Liposome

Form

Liquid

Purification

None

Buffer

25 mM Tris-HCl of pH8.0 containing 2% glycerol

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

FCAMR

Full Name

Fc fragment of IgA and IgM receptor

Introduction

FCAMR (Fc Fragment Of IgA And IgM Receptor) is a Protein Coding gene. Among its related pathways are Cell surface interactions at the vascular wall and Response to elevated platelet cytosolic Ca²⁺. Gene Ontology (GO) annotations related to this gene include transmembrane signaling receptor activity and IgM binding. An important paralog of this gene is PIGR

Alternative Names

FCA/MR; FKSG87; Fc alpha/mu receptor,immunity related factor,receptor for Fc fragment of IgA and IgM

Gene ID

[83953](#)

UniProt ID

[Q8WWV6](#)