

Product Information

MemDX™ Membrane Protein Human FMR1NB (FMR1 neighbor) for Antibody Discovery

Cat. No.: **MP0371X**

This product is for research use only and is not intended for diagnostic use.

This product is a 53.57 kDa Human FMR1NB membrane protein expressed in *in vitro* wheat germ expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

FMR1NB

Protein Length

Full-length

Molecular Weight

53.57 kDa

TMD

2

Sequence

MSSHRRKAKGRNRRSHRAMRVAHLELATYELAATESNPESHHPGYEAAMADRPQPGWRESLKMRVSKPFGMLMLSIWILLFVCY

Product Description

Application

Enzyme-linked Immunoabsorbent Assay, Western Blot (Recombinant protein), Antibody Production, Protein Array

Expression Systems

in vitro wheat germ expression system

Tag

GST-tag at N-terminal

Form

Liquid

Purification

Glutathione Sepharose 4 Fast Flow

Buffer

50 mM Tris-HCl, 10 mM reduced Glutathione, pH=8.0 in the elution buffer

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

FMR1NB

Full Name

FMR1 neighbor

Introduction

FMR1NB (FMR1 Neighbor) is a Protein Coding gene. Diseases associated with FMR1NB include Splenic Artery Aneurysm and Sarcoma

Alternative Names

FLJ25736; NY-SAR-35; OTTHUMP00000024202

Gene ID

[158521](#)

UniProt ID

[Q8N0W7](#)