

Product Information

MemDX™ Antibody Discovery - Human SLAMF7 / CRACC / CD319 (23-226) Membrane

Protein, Partial, -His tag, [FITC]

Cat. No.: **MP0410F**

This product is for research use only and is not intended for diagnostic use.

This membrane protein is Human SLAMF7 / CRACC / CD319 (23-226). It has been tested in SDS-PAGE, FACS. We provide this protein to facilitate your membrane protein antibody discovery and development.

Product Specifications

Host Species

Human

Target Protein

SLAMF7 / CRACC / CD319

Protein Length

ECD

Molecular Weight

The protein has a calculated MW of 23.2 kDa. The protein migrates as 33-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Sequence

AA Ser 23 - Met 226 (Accession # AAH27867).

Product Description

Activity

Yes

Application

SDS-PAGE, FACS

Expression Systems

HEK293

Tag

His tag at the C-terminus

Protein Format

Soluble

Form

LYOPH

Reconstitution

Please see Certificate of Analysis for specific instructions.

Endotoxin

<1.0 EU/μg by the LAL method

Conjugation

FITC

Purity

>90% as determined by SDS-PAGE.

Buffer

Lyophilized from 0.22 μm filtered solution in PBS, pH7.4. Normally trehalose is added as protectant before lyophilization.

Storage

Stored at lyophilized form at -20°C or lower. Avoid repeated freeze-thaw cycles.

The antigen can be stable for 12 months in lyophilized form after storage at -20°C to -80°C, 3 months under sterile conditions after reconstitution after storage at -80°C.

Target

Target Protein

SLAMF7 / CRACC / CD319

Full Name

SLAM family member 7

Alternative Names

19A; CS1; CD319; CRACC; SLAM family member 7; 19A24 protein; CD2 subset 1; CD2-like receptor activating cytotoxic cells; membrane protein FOAP-12; novel LY9 (lymphocyte antigen 9) like protein; protein 19A

Gene ID

[57823](#)

UniProt ID

[Q9NQ25](#)