

Product Information

MemDX™ Membrane Protein Human GPR84 (G protein-coupled receptor 84) for Antibody

Discovery

Cat. No.: **MP0475X**

This product is for research use only and is not intended for diagnostic use.

This product is a 43.7 kDa Human GPR84 membrane protein expressed in *in vitro* wheat germ expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

GPR84

Protein Length

Full-length

Molecular Weight

43.7 kDa

TMD

7

Sequence

MWNSSDANFSCYHESVLGYRYVAVSWGVVAVTGTGTVGNVLTLLALAIQPKLRTRFNLLIANLTLADLLYCTLLQPFSVDTYLHLHWR

Product Description

Application

Antibody Production

Expression Systems

in vitro wheat germ expression system

Tag

NO

Protein Format

Liposome

Form

Liquid

Purification

None

Buffer

25 mM Tris-HCl of pH8.0 containing 2% glycerol

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

GPR84

Full Name

G protein-coupled receptor 84

Introduction

The GPR84 gene is conserved in chimpanzee, Rhesus monkey, dog, cow, mouse, rat, zebrafish, fruit fly, mosquito, and frog.

Alternative Names

EX33; GPCR4

Gene ID

[53831](#)

UniProt ID

[Q9NQS5](#)