

Product Information

MemDX™ Antibody Discovery - Cynomolgus / Rhesus macaque TIGIT (89-208) Membrane Protein, Partial, -hlgG1 Fc tag

Cat. No.: **MP0483F**

This product is for research use only and is not intended for diagnostic use.

This membrane protein is Cynomolgus / Rhesus macaque TIGIT (89-208). It has been tested in SDS-PAGE, ELISA, BLI. We provide this protein to facilitate your membrane protein antibody discovery and development.

Product Specifications

Host Species

Cynomolgus / Rhesus macaque

Target Protein

TIGIT

Protein Length

ECD

Molecular Weight

The protein has a calculated MW of 39.9 kDa. The protein migrates as 45-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Sequence

AA Met 89 - Ile 208 (Accession # G7NXM4-1).

Product Description

Activity

Yes

Application

SDS-PAGE, ELISA, BLI

Expression Systems

HEK293

Tag

Human IgG1 Fc tag at the C-terminus

Protein Format

Soluble

Form

LYOPH

Reconstitution

Please see Certificate of Analysis for specific instructions.

Endotoxin

<1.0 EU/μg by the LAL method

Purity

>95% as determined by reduced SDS-PAGE.

Buffer

Lyophilized from 0.22 μm filtered solution in 50 mM Tris, 100 mM Glycine, pH7.5. Normally trehalose is added as protectant before lyophilization.

Storage

Stored at lyophilized form at -20°C or lower. Avoid repeated freeze-thaw cycles.

The antigen can be stable for 12 months in lyophilized form after storage at -20°C to -80°C, 3 months under sterile conditions after reconstitution after storage at -80°C.

Target

Target Protein

TIGIT

Full Name

T cell immunoreceptor with Ig and ITIM domains

Introduction

This gene encodes a member of the PVR (poliovirus receptor) family of immunoglobulin proteins. The product of this gene is expressed on several classes of T cells including follicular B helper T cells (TFH). The protein has been shown to bind PVR with high affinity; this binding is thought to assist interactions between TFH and dendritic cells to regulate T cell dependent B cell responses

Alternative Names

TIGIT, VSIG9, VSTM3

Gene ID

[102121588](#)

UniProt ID

[A0A2K5UW92](#)