

Product Information

MemDX™ Membrane Protein Human LRRN4CL (LRRN4 C-terminal like) for Antibody

Discovery

Cat. No.: **MP0644X**

This product is for research use only and is not intended for diagnostic use.

This product is a 51.7 kDa Human LRRN4CL membrane protein expressed in *in vitro* wheat germ expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

LRRN4CL

Protein Length

Full-length

Molecular Weight

51.7 kDa

TMD

1

Sequence

MLGSPCLLWLLAVTFLVPRAQLAPQDFEEEEADETETAWPPLPAVPCDYDHCRLQVPCKELQRVGPAACLCPLGSSPAQPPDF

Product Description

Application

Enzyme-linked Immunoabsorbent Assay, Western Blot (Recombinant protein), Antibody Production, Protein Array

Expression Systems

in vitro wheat germ expression system

Tag

GST-tag at N-terminal

Form

Liquid

Purification

Glutathione Sepharose 4 Fast Flow

Buffer

50 mM Tris-HCl, 10 mM reduced Glutathione, pH=8.0 in the elution buffer

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

LRRN4CL

Full Name

LRRN4 C-terminal like

Introduction

The LRRN4CL gene is conserved in chimpanzee, Rhesus monkey, dog, cow, mouse, rat, and frog.

Alternative Names

MGC61707

Gene ID

[221091](#)

UniProt ID

[Q8ND94](#)