

# Product Information

## MemDX™ Membrane Protein Human SLC2A4 (Solute carrier family 2 member 4) for Antibody Discovery

Cat. No.: **MP0646J**

This product is for research use only and is not intended for diagnostic use.

This product is a 54.6 kDa Human SLC2A4 membrane protein expressed in HEK293T. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

SLC2A4

#### Protein Length

Full-length

#### Protein Class

Druggable Genome, Transmembrane

#### Molecular Weight

54.6 kDa

#### TMD

12

#### Sequence

MPSGFQQIGSEDPQQQRVTGTLVLAVFSAVLGSLQFGYNIGVINAPQKVIEQSYNETWLGRQGPEGPS  
SIPPGTLTTLWALSVAIFSVGGMISSFLIGISQWLGRKRAMLVNNVLAVLGGSLMGLANAAASYEMLIL  
GRFLIGAYSGLTSGLVPMYVGEIAPTHLRGALGTNLQLAIVIGILIAQVLGLESLLGTASLWPLLLGLTV  
LPALLQLVLLPFCPESPRYLYIIQNLEGPARKSLKRLTGWADVSGVLAELKDEKRLERERPLSLLQLLG  
SRTHRQPLIIAVVLQLSQQLSGINAVFYYSTSIFETAGVGQPAYATIGAGVVNTVFTLVSVLLVERAGRR  
TLHLLGLAGMCGCAILMTVALLLLERVPAISYVSIVAIFGFVAFGEIGPGPIPWFIWAELFSQGPRPAAM  
AVAGFSNWTNFIIGMGFQYVAEAMGPYVFLFAVLLLGFFIFTFLRVPETRGRTFDQISAAFHRTPSLL  
EQEVKPSTELEYLGPDEND

### Product Description

#### Expression Systems

HEK293T

#### Tag

C-Myc/DDK

**Form**

Liquid

**Purification**

Anti-DDK affinity column followed by conventional chromatography steps

**Purity**

> 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer**

25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

**Storage**

Store at +4°C for up to one week or several months at -80°C

**Target****Target Protein**

SLC2A4

**Full Name**

Solute carrier family 2 member 4

**Introduction**

This gene is a member of the solute carrier family 2 (facilitated glucose transporter) family and encodes a protein that functions as an insulin-regulated facilitative glucose transporter. In the absence of insulin, this integral membrane protein is sequestered within the cells of muscle and adipose tissue. Within minutes of insulin stimulation, the protein moves to the cell surface and begins to transport glucose across the cell membrane. Mutations in this gene have been associated with noninsulin-dependent diabetes mellitus (NIDDM).

**Alternative Names**

GLUT4

**Gene ID**

[6517](#)

**UniProt ID**

[P14672](#)