

Product Information

MemDX™ Membrane Protein Human IRAG2 (Inositol 1,4,5-triphosphate receptor associated 2) for Antibody Discovery

Cat. No.: **MP0652X**

This product is for research use only and is not intended for diagnostic use.

This product is a 54.9 kDa Human IRAG2 membrane protein expressed in *in vitro* wheat germ expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

IRAG2

Protein Length

Full-length

Molecular Weight

54.9 kDa

TMD

1

Sequence

MNDDPSMEENGVERVCPESELLQSREYSSLPLPRHTSSTDGTITSSDPGLEILNMASCDLDRNSLCKKEEDTRSASPTIEAQGTSPAA

Product Description

Application

Antibody Production

Expression Systems

in vitro wheat germ expression system

Tag

NO

Protein Format

Liposome

Form

Liquid

Purification

None

Buffer

25 mM Tris-HCl of pH8.0 containing 2% glycerol

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

IRAG2

Full Name

Inositol 1,4,5-triphosphate receptor associated 2

Introduction

The protein encoded by this gene is expressed in a developmentally regulated manner in lymphoid cell lines and tissues. The protein is localized to the cytoplasmic face of the endoplasmic reticulum.

Alternative Names

JAW1; LRMP

Gene ID

[4033](#)

UniProt ID

[Q12912](#)