

Product Information

MemDX™ Antibody Discovery - Human Nectin-2 / CD112 (32-360) Membrane Protein, Partial, -hIgG1 Fc -Avi tag, [Biotin]

Cat. No.: **MP0693F**

This product is for research use only and is not intended for diagnostic use.

This membrane protein is Human Nectin-2 / CD112 (32-360). It has been tested in SDS-PAGE, ELISA. We provide this protein to facilitate your membrane protein antibody discovery and development.

Product Specifications

Host Species

Human

Target Protein

Nectin-2 / CD112

Protein Length

ECD

Molecular Weight

The protein has a calculated MW of 64.0 kDa. As a result of glycosylation, the protein migrates as 66 kDa under reducing (R) condition, and 130 kDa under non-reducing (NR) condition (SDS-PAGE).

Sequence

AA Gln 32 - Leu 360 (Accession # Q92692-2).

Product Description

Activity

Yes

Application

SDS-PAGE, ELISA

Expression Systems

HEK293

Tag

Human IgG1 Fc tag at the C-terminus, followed by a Avi tag

Protein Format

Soluble

Form

LYOPH

Reconstitution

Please see Certificate of Analysis for specific instructions.

Endotoxin

<1.0 EU/μg by the LAL method

Conjugation

Biotin

Purity

>95% as determined by SDS-PAGE.

Buffer

Lyophilized from 0.22 μm filtered solution in Tris with Glycine, Arginine and NaCl, pH7.5. Normally trehalose is added as protectant before lyophilization.

Storage

Stored at lyophilized form at -20°C or lower. Avoid repeated freeze-thaw cycles.

The antigen can be stable for 12 months in lyophilized form after storage at -20°C to -80°C, 3 months under sterile conditions after reconstitution after storage at -80°C.

Target

Target Protein

Nectin-2 / CD112

Full Name

nectin cell adhesion molecule 2

Introduction

This gene encodes a single-pass type I membrane glycoprotein with two Ig-like C2-type domains and an Ig-like V-type domain. This protein is one of the plasma membrane components of adherens junctions. It also serves as an entry for certain mutant strains of herpes simplex virus and pseudorabies virus, and it is involved in cell to cell spreading of these viruses. Variations in this gene have been associated with differences in the severity of multiple sclerosis. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.

Alternative Names

HVEB; PRR2; CD112; PVRL2; PVRR2; nectin-2; herpesvirus entry protein B; poliovirus receptor-like 2; poliovirus receptor-related 2 (herpesvirus entry mediator B)

Gene ID

[5819](#)

UniProt ID

[Q92692](#)