

Product Information

MemDX™ Membrane Protein Human PRLHR (Prolactin releasing hormone receptor) without tag for Antibody Discovery

Cat. No.: **MP1043X**

This product is for research use only and is not intended for diagnostic use.

This product is a 41.1 kDa Human PRLHR membrane protein expressed in *in vitro* wheat germ expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

PRLHR

Protein Length

Full-length

Molecular Weight

41.1 kDa

TMD

7

Sequence

MASSTTRGPRVSDLFSGLPVAVTTPANQSAEASAGNGSVAGADAPAVTPFQSLQLVHQLKGLIVLLYSVVVVGLVGNCLLVLIAR

Product Description

Application

Antibody Production

Expression Systems

in vitro wheat germ expression system

Tag

NO

Protein Format

Liposome

Form

Liquid

Buffer

25 mM Tris-HCl of pH8.0 containing 2% glycerol

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

PRLHR

Full Name

Prolactin releasing hormone receptor

Introduction

PRLHR is a 7-transmembrane domain receptor for prolactin-releasing hormone (PRLH; MIM 602663) that is highly expressed in anterior pituitary.

Alternative Names

GR3; GPR10; PrRPR; prolactin-releasing peptide receptor; G-protein coupled receptor 10; hGR3; prRP receptor

Gene ID

[2834](#)

UniProt ID

[P49683](#)