

Product Information

MemDX™ Antibody Discovery - Human BTLA (31-150) (31-150) Membrane Protein, Partial, - hIgG1 Fc tag

Cat. No.: **MP1057F**

This product is for research use only and is not intended for diagnostic use.

This membrane protein is Human BTLA (31-150) (31-150). It has been tested in SDS-PAGE, ELISA. We provide this protein to facilitate your membrane protein antibody discovery and development.

Product Specifications

Host Species

Human

Target Protein

BTLA (31-150)

Protein Length

ECD

Molecular Weight

The protein has a calculated MW of 40.4 kDa. The protein migrates as 46-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Sequence

AA Lys 31 - Ser 150 (Accession # AAP44003.1).

Product Description

Activity

Yes

Application

SDS-PAGE, ELISA

Expression Systems

HEK293

Tag

Human IgG1 Fc tag at the C-terminus

Protein Format

Soluble

Form

LYOPH

Reconstitution

Please see Certificate of Analysis for specific instructions.

Endotoxin

<1.0 EU/μg by the LAL method

Purity

>92% as determined by SDS-PAGE.

Buffer

Lyophilized from 0.22 μm filtered solution in Tris with Glycine, Arginine and NaCl, pH7.5. Normally trehalose is added as protectant before lyophilization.

Storage

Stored at lyophilized form at -20°C or lower. Avoid repeated freeze-thaw cycles.

The antigen can be stable for 12 months in lyophilized form after storage at -20°C to -80°C, 3 months under sterile conditions after reconstitution after storage at -80°C.

Target

Target Protein

BTLA (31-150)

Full Name

B and T lymphocyte associated

Introduction

This gene encodes a member of the immunoglobulin superfamily. The encoded protein contains a single immunoglobulin (Ig) domain and is a receptor that relays inhibitory signals to suppress the immune response. Alternative splicing results in multiple transcript variants. Polymorphisms in this gene have been associated with an increased risk of rheumatoid arthritis.

Alternative Names

BTLA1; CD272; B- and T-lymphocyte attenuator; B- and T-lymphocyte-associated protein

Gene ID

[151888](#)

UniProt ID

[Q7Z6A9](#)