

Product Information

MemDX™ Antibody Discovery - Rhesus macaque CD155 / PVR (28-343) Membrane Protein, Partial, -His tag

Cat. No.: **MP1070F**

This product is for research use only and is not intended for diagnostic use.

This membrane protein is Rhesus macaque CD155 / PVR (28-343). It has been tested in SDS-PAGE, ELISA, BLI. We provide this protein to facilitate your membrane protein antibody discovery and development.

Product Specifications

Host Species

Rhesus macaque

Target Protein

CD155 / PVR

Protein Length

ECD

Molecular Weight

The protein has a calculated MW of 36.3 kDa. The protein migrates as 48-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Sequence

AA Asp 28 - Asn 343 (Accession # Q0MSE6-1).

Product Description

Activity

Yes

Application

SDS-PAGE, ELISA, BLI

Expression Systems

HEK293

Tag

His tag at the C-terminus

Protein Format

Soluble

Form

LYOPH

Reconstitution

Please see Certificate of Analysis for specific instructions.

Endotoxin

<1.0 EU/μg by the LAL method

Purity

>90% as determined by SDS-PAGE.

Buffer

Lyophilized from 0.22 μm filtered solution in PBS, pH7.4. Normally trehalose is added as protectant before lyophilization.

Storage

Stored at lyophilized form at -20°C or lower. Avoid repeated freeze-thaw cycles.

The antigen can be stable for 12 months in lyophilized form after storage at -20°C to -80°C, 3 months under sterile conditions after reconstitution after storage at -80°C.

Target

Target Protein

CD155 / PVR

Full Name

poliovirus receptor

Introduction

The protein encoded by this gene is a transmembrane glycoprotein belonging to the immunoglobulin superfamily. The external domain mediates cell attachment to the extracellular matrix molecule vitronectin, while its intracellular domain interacts with the dynein light chain Tctex-1/DYNLT1. The gene is specific to the primate lineage, and serves as a cellular receptor for poliovirus in the first step of poliovirus replication. Multiple transcript variants encoding different isoforms have been found for this gene.

Alternative Names

poliovirus receptor, nectin-like protein 5

Gene ID

[712851](#)

UniProt ID

[Q0MSE6](#)