

# Product Information

## **MemDX™ Membrane Protein Human SEZ6L (Seizure related 6 homolog like) expressed in *In vitro* wheat germ expression system for Antibody Discovery**

Cat. No.: **MP1122X**

This product is for research use only and is not intended for diagnostic use.

This product is a 130.68 kDa Human SEZ6L membrane protein expressed in *In vitro* wheat germ expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

SEZ6L

#### Protein Length

Full-length

#### Molecular Weight

130.68 kDa

#### TMD

1

#### Sequence

MPAARPPAAGLRGISLFLALLLGSPAAALERDALPEGDASPLGPYLLPSGAPERGSPGKEHPEERVVTAPPSSSQSAEVLGELVLD

### Product Description

#### Application

Enzyme-linked Immunoabsorbent Assay, Western Blot (Recombinant protein), Antibody Production, Protein Array

#### Expression Systems

*In vitro* wheat germ expression system

#### Tag

GST-tag at N-terminal

#### Form

Liquid

#### Purification

Glutathione Sepharose 4 Fast Flow

**Buffer**

50 mM Tris-HCl, 10 mM reduced Glutathione, pH=8.0

**Storage**

Store at +4°C for up to one week or several months at -80°C

**Target****Target Protein**

SEZ6L

**Full Name**

Seizure related 6 homolog like

**Introduction**

May contribute to specialized endoplasmic reticulum functions in neurons.

**Alternative Names**

seizure 6-like protein; seizure related gene 6-like

**Gene ID**

[23544](#)

**UniProt ID**

[Q9BYH1](#)