

Product Information

MemDX™ Membrane Protein Human CD44 (CD44 molecule (Indian blood group)) for

Antibody Discovery

Cat. No.: **MP1152J**

This product is for research use only and is not intended for diagnostic use.

This product is a 79.2 kDa Human CD44 membrane protein expressed in HEK293T. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

CD44

Protein Length

Full-length

Protein Class

Adult stem cells, Cancer stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Stem cell relevant signaling - DSL/Notch pathway, Transmembrane

Molecular Weight

79.2 kDa

TMD

1

Sequence

MDKFWWWHAAGWGLCLVPLSLAQIDLNITCRFAGVFHVEKNGRYSISRTEAADLCKAFNSTLPTMAQMEKAL
SIGFETCRYGFIIEGHVVIPRIHPNSICAANNTGVYILTSNTSQYDTYCFNASAPPEEDCTSVTDLPNAFD
GPITITIVNRDGTTRYVQKGEYRTNPEDIYPSNPTDDDVSSGSSSERSSTSGGYIFYTFSTVHPIPEDDSP
WITDSTDRIPTATLMSTSATATETATKRQETWDWFSWFLPSESKNHLHTTTQMAGTSSNTISAGWEPNE
ENEDERDRHLSFSGSGIDDDFISSTISTTPRAFDHTKQNQDWTQWNPSHSNPEVLLQTTTRMTDVRN
GTTAYEGNWNPEAHPLIHHEHHEEEETPHSTSTIQATPSSTTEETATQKEQWFGNRWHEGYRQTPKEDS
HSTTGTAASAHTSHPMQGRTPSPEDSSWTDFFNPISHPMGRGHQAGRRMDMDSSHSITLQPTANPNTG
LVEDLDRTGPLSMTTQQSNSQSFSTSHGLEEDKDHPPTSTLTSSNRNDVTGGRRDPNHSEGSTTLLEGY
TSHYPHTKESRTFIPVTSAKTGSFGVTAVTVGDSNSNVNRSLSGDQDTFHPSGGSHHTHGSESDGHSHG
QEGGANTTSGPIRTPQIPEWLILASLLALILAVCIAVNSRRRCGQKKLVINSGNGAVEDRKPSGLN
GEASKSQEMVHLVNKESSETPDQFMTADETRNLQNVDKIGV

Product Description

Expression Systems

HEK293T

Tag

C-Myc/DDK

Form

Liquid

Purification

Anti-DDK affinity column followed by conventional chromatography steps

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer

25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

CD44

Full Name

CD44 molecule (Indian blood group)

Introduction

The protein encoded by this gene is a cell-surface glycoprotein involved in cell-cell interactions, cell adhesion and migration. It is a receptor for hyaluronic acid (HA) and can also interact with other ligands, such as osteopontin, collagens, and matrix metalloproteinases (MMPs). This protein participates in a wide variety of cellular functions including lymphocyte activation, recirculation and homing, hematopoiesis, and tumor metastasis. Transcripts for this gene undergo complex alternative splicing that results in many functionally distinct isoforms, however, the full length nature of some of these variants has not been determined. Alternative splicing is the basis for the structural and functional diversity of this protein, and may be related to tumor metastasis.

Alternative Names

IN; LHR; MC56; MDU2; MDU3; MIC4; Pgp1; CDW44; CSPG8; HCELL; HUTCH-I; ECMR-III

Gene ID

[960](#)

UniProt ID

[P16070](#)