

Product Information

MemDX™ Antibody Discovery - Human CD31 / PECAM-1 (28-601) Membrane Protein, Partial, -His tag, [PE]

Cat. No.: **MP1184F**

This product is for research use only and is not intended for diagnostic use.

This membrane protein is Human CD31 / PECAM-1 (28-601). It has been tested in SDS-PAGE. We provide this protein to facilitate your membrane protein antibody discovery and development.

Product Specifications

Host Species

Human

Target Protein

CD31 / PECAM-1

Protein Length

ECD

Molecular Weight

The protein has a calculated MW of 66.4 kDa. The protein migrates as 80-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Sequence

AA Gln 28 - Lys 601 (Accession # P16284-1).

Product Description

Application

SDS-PAGE

Expression Systems

HEK293

Tag

His tag at the C-terminus

Protein Format

Soluble

Form

LYOPH

Reconstitution

Please see Certificate of Analysis for specific instructions.

Endotoxin

<1.0 EU/μg by the LAL method

Conjugation

PE

Purity

>95% as determined by SDS-PAGE.

Buffer

Lyophilized from 0.22 μm filtered solution in 50 mM Tris, 150 mM NaCl, pH7.5 . Normally trehalose is added as protectant before lyophilization.

Storage

Stored at lyophilized form at -20°C or lower. Avoid repeated freeze-thaw cycles.

The antigen can be stable for 12 months in lyophilized form after storage at -20°C to -80°C, 3 months under sterile conditions after reconstitution after storage at -80°C.

Target**Target Protein**

CD31 / PECAM-1

Full Name

platelet and endothelial cell adhesion molecule 1

Introduction

The protein encoded by this gene is found on the surface of platelets, monocytes, neutrophils, and some types of T-cells, and makes up a large portion of endothelial cell intercellular junctions. The encoded protein is a member of the immunoglobulin superfamily and is likely involved in leukocyte migration, angiogenesis, and integrin activation.

Alternative Names

CD31; PECA1; GPIIA'; PECAM-1; endoCAM; CD31/EndoCAM; platelet endothelial cell adhesion molecule; CD31 antigen; platelet endothelial cell adhesion molecule-1

Gene ID

[5175](#)

UniProt ID

[P16284](#)