

Product Information

MemDX™ Membrane Protein Human ERBB2 (Erb-b2 receptor tyrosine kinase 2) expressed in sf9 cells for Antibody Discovery

Cat. No.: **MP1309J**

This product is for research use only and is not intended for diagnostic use.

This product is a 64 kDa Human ERBB2 membrane protein expressed in Sf9. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ERBB2

Protein Length

Partial (676-1255aa)

Protein Class

Druggable Genome, Protein Kinase, Transmembrane

Molecular Weight

64 kDa

TMD

1

Sequence

MELAALCRWGLLLALLPPGAASTQVCTGTMKLRLPASPETHLDMLRHL YQGCQVVQGNLELT YLPTNAS
LSFLQDIQEVQGYVLI AHNQVRQVPLQRLRIVRG TQLFEDNYALAVLDNGDPLNNTTPVTGASPGGLREL
QLRSLTEILKGGVLIQRNPQLCYQDTILWKDIFHKNNQLAL TLIDTNRSRACHPCSPMCKGSR CWGESSE
DCQSLTRTV CAGGCARCKGPLPTDCCHEQCAAGCTGPKHSDCLACLHFNHSGICELHCPALVTYNTDTFE
SMPNPEG RYTFGASCVTACPYNYLSTDVGSC TLVCP LHNQEVTAEDGTQRCEKCSKPCARVCYGLGMEHL
REVRAVTSANIQEFAGCKKIFGSLAFLPESFDGDPASNTAPLQPEQLQVFETLEEITGYLYISAWPDSL P
DLSVFQNLQVIRGRILHNGAYSLTLQGLGISWLGLRSLRELGSG LALIHNNTHLCFVHTVPWDQLFRNPH
QALLHTANRPEDEC VGEGLACHQLCARGHCWGP GPTQCVNCSQFLRGQECVEECRVLQGLPREYVNARHC
LPCHPEQC PQNGSVTCFGL EADQC VACA HYKDP PFCVARCP SGVKPDL SYMPIWKFPDEEGACQPCPINC
THSCVDLDDKCPAEQRASPLTSIISAVVGILLVVVLGVVFGILIKRRQQKIRKYTMRRLLQETELVEPL
TPSGAMPNQAQMRILKETELRKVKVLGSGAFGT VYKGIWIPDGENVKIPVAIKVLRENTSPKANKEILDE
AYVMAGVGSPYVSRLLGICLTSTVQLVTQLMPYGC LLDHVREN RGR LGSQDLLNWC MQIAKGMSYLEDVR
LVHRDLAARNVLVKSPNHVKITDFGLARLLDIDET EYHADGGKVPIKWMAL ESILRRRFTHQSDVWSYGV
TVWELMTFGAKPYDGIPAREIPDLLEKGERLPQPPICTIDVYMIMVKCWMIDSECRPRFREL VSEFSRMA
RDPQRFVVIQNE DLGPASPLDSTFYRSLLEDDDMGDLVDAEEYLVPQQGFFCPDPAPGAGGMVHHRHRS
STRSGGGDLTLGLEPSEEEAPRSP LAPSEGAGSDVFDGDLGMGA AKGLQSLPTHDPSP LQRYSEDPTVPL
PSETDGYVAPLTCSPQPEYVNQPDVRPQPPSPREGPLPAARPAGATLERPKTLSPGKNGVVKDVFAFGGA
VENPEYLT PQGGAAPQPHPPPAFSPA FDNLYYWDQDPPERGAPPSTFKGTPTAENPEYLGLDVPV

Product Description

Expression Systems

Sf9

Tag

C-DDK

Form

Liquid

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer

50mM Tris-HCl pH8.0, 150mM NaCl, 20% glycerol

Storage

Store at +4°C for up to one week or several months at -80°C

Target

Target Protein

ERBB2

Full Name

Erb-b2 receptor tyrosine kinase 2

Introduction

This gene encodes a member of the epidermal growth factor (EGF) receptor family of receptor tyrosine kinases. This protein has no ligand binding domain of its own and therefore cannot bind growth factors. However, it does bind tightly to other ligand-bound EGF receptor family members to form a heterodimer, stabilizing ligand binding and enhancing kinase-mediated activation of downstream signalling pathways, such as those involving mitogen-activated protein kinase and phosphatidylinositol-3 kinase. Allelic variations at amino acid positions 654 and 655 of isoform a (positions 624 and 625 of isoform b) have been reported, with the most common allele, Ile654/Ile655, shown here. Amplification and/or overexpression of this gene has been reported in numerous cancers, including breast and ovarian tumors. Alternative splicing results in several additional transcript variants, some encoding different isoforms and others that have not been fully characterized.

Alternative Names

CD340; HER-2; HER-2/neu; HER2; MLN 19; NEU; NGL; TKR1

Gene ID

[2064](#)

UniProt ID

[P04626](#)