

Product Information

MemDX™ Membrane Protein Human TICAM2 (Toll like receptor adaptor molecule 2)

expressed in *In vitro* wheat germ expression system for Antibody Discovery

Cat. No.: **MP1356X**

This product is for research use only and is not intended for diagnostic use.

This product is a 53.3 kDa Human TICAM2 membrane protein expressed in *In vitro* wheat germ expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TICAM2

Protein Length

Full-length

Molecular Weight

53.3 kDa

TMD

1

Sequence

MGIGKSKINSCPLSLSWGKRHSVDTSPGYHESDSKKSEDLSLCNVAEHSNTTEGPTGKQEGAQSVEEMFEEEEAEVFLKFVILHA

Product Description

Application

Enzyme-linked Immunoabsorbent Assay, Western Blot (Recombinant protein), Antibody Production, Protein Array

Expression Systems

in vitro wheat germ expression system

Tag

GST-tag at N-terminal

Protein Format

Liposome

Form

Liquid

Purification

Glutathione Sepharose 4 Fast Flow

Buffer

50 mM Tris-HCl, 10 mM reduced Glutathione, pH=8.0

Storage

Store at +4°C for up to one week or several months at -80°C

Target

Target Protein

TICAM2

Full Name

Toll like receptor adaptor molecule 2

Introduction

TIRP is a Toll/interleukin-1 receptor (IL1R; MIM 147810) (TIR) domain-containing adaptor protein involved in Toll receptor signaling.

Alternative Names

TIRP; TRAM; TIRAP3; MyD88-4; TICAM-2; TIR domain-containing adapter molecule 2; NF-kappa-B-activating protein 502; TRIF-related adaptor molecule; cytoplasmic adaptor; putative NF-kappa-B-activating protein 502; toll-like receptor adaptor protein 3; toll/interleukin-1 receptor (TIR) domain-containing adapter protein

Gene ID

[353376](#)

UniProt ID

[Q86XR7](#)