

Product Information

MemDX™ Membrane Protein Human RYR3 (Ryanodine receptor 3) for Antibody Discovery

Cat. No.: **MP1406J**

This product is for research use only and is not intended for diagnostic use.

This product is a 32.5 kDa Human RYR3 membrane protein expressed in E.coli. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

RYR3

Protein Length

Partial (3934-4181aa)

Protein Class

Ion Channel

Molecular Weight

32.5 kDa

Sequence

DGKGIISKKEFQKAMEGQKQYTSQSEIDFLLSCAEADENDMFNYVDFVDRFHPEAKDIGFNVAVLLTNLSEHMPNDSRLKCLLDPAES

Product Description

Expression Systems

E.coli

Tag

N-6xHis

Form

Liquid or Lyophilized powder

Reconstitution

Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration).

Purity

>90% as determined by SDS-PAGE

Buffer

Liquid: Tris/PBS-based buffer, 5%-50% glycerol

Lyophilized powder: Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

RYR3

Full Name

Ryanodine receptor 3

Introduction

The protein encoded by this gene is a ryanodine receptor, which functions to release calcium from intracellular storage for use in many cellular processes. For example, the encoded protein is involved in skeletal muscle contraction by releasing calcium from the sarcoplasmic reticulum followed by depolarization of T-tubules. Two transcript variants encoding different isoforms have been found for this gene.

Alternative Names

Brain ryanodine receptor calcium release channel; Brain ryanodine receptor-calcium release channel; Brain type ryanodine receptor; Brain-type ryanodine receptor; HBRR; Ryanodine receptor 3; RYR 3; RYR-3; RyR3; RYR3_HUMAN

Gene ID

[6263](#)

UniProt ID

[Q15413](#)