

Product Information

MemDX™ Membrane Protein Human FPR3 (Formyl peptide receptor 3) Full Length

Cat. No.: **MPC0095K**

This product is for research use only and is not intended for diagnostic use.

This product is a 39.9 kDa Human FPR3 membrane protein expressed in Baculovirus/Insect expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

FPR3

Protein Length

Full length

Protein Class

GPCR

Molecular Weight

39.9 kDa

TMD

7

Sequence

METNFSIPLNETEEVLPEPAGHTVLWIFSLLVHGVTFVFGVLGNGLVIWV
AGFRMTRTVNTICYLNLALADFSFSAILPFRMVSVAMREKWPFGSFLCKL
VHVMIDINLFVSVYLITIIALDRCICVLHPAWAQNHRTMSLAKRVMTGLW
IFTIVLTLPNFIWTTISTTNGDTYCFNFAFWGDTAVERLNVFITMAKV
FLILHFIIGFSVPMSIITVCYGIIAAKIHRNHMIKSSRPLRVFAAVVASF
FICWFPYELIGILMAVWLKEMLLNGKYKIILVLINPTSSLAFFNSCLNPI
LYVFMGRNFQERLIRSLPTSLERALTEVPDSAQTSNTDTTSASPPEETEL
QAM

Product Description

Expression Systems

Baculovirus/Insect expression system

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

FPR3

Full Name

Formyl peptide receptor 3

Alternative Names

FMLPY; FPRH1; FPRH2; FPRL2; RMLP-R-I; FMLP-R-II; FML2_HUMAN; N-formyl peptide receptor 3; FMLP-related receptor II; formyl peptide receptor-like 2; FPR3; Formyl peptide receptor 3

Gene ID

[2359](#)

UniProt ID

[P25089](#)