

Product Information

MemDX™ Membrane Protein Human GPR152 (G protein-coupled receptor 152) Full Length

Cat. No.: **MPC0130K**

This product is for research use only and is not intended for diagnostic use.

This product is a 50.9 kDa Human GPR152 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

GPR152

Protein Length

Full length

Protein Class

GPCR

Molecular Weight

50.9 kDa

TMD

7

Sequence

MDTTMEADLGATGHRPRTLEDDEDSYPQGGWDTVFLVALLLLGLPANGLM
AWLAGSQARHGAGTRLALLLLSLALSDFLFLAAAFQILEIRHGGHWPLG
TAACRFYYFLWGVSYSSGLFLLAALSLDRCLLALCPHWYPGHRPVRLPLW
VCAGVWVLATLFSVPWLVFPEAAVWWYDLVICLDFWDSEELSLRMLEVLG
GFLPFLLLLCHVLTQATACRTCHRQQQPAACRGFARVARTILSAYVVLR
LPYQLAQLLYLAFLWDVYSGYLLWEALVYSYDYLILLNSCLSPFLCLMASA
DLRTLRLSVLSSFAAALCEERPGSFTPTPTQQLDSEGPTLPEPMAEAQS
QMDPVAAQPQVNPTLQPRSDPTAQPQLNPTAQPQSDPTAQPQLNLMAQPQS
DSVAQPQADTNVQTPAPAASSVPSPCDEASPTPSSHPTPGALEDPATPPA
SEGESPSSTPPEAAPGAGPT

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

GPR152

Full Name

G protein-coupled receptor 152

Alternative Names

PGR5; probable G-protein coupled receptor 152; G protein-coupled receptor PGR5; GPR152; G protein-coupled receptor 152

Gene ID

[390212](#)

UniProt ID

[Q8TDT2](#)