

Product Information

MemDX™ Membrane Protein Human GPR17 (G protein-coupled receptor 17) Full Length

Cat. No.: **MPC0138K**

This product is for research use only and is not intended for diagnostic use.

This product is a 40.9 kDa Human GPR17 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

GPR17

Protein Length

Full length

Protein Class

GPCR

Molecular Weight

40.9 kDa

TMD

7

Sequence

MSKRSWWAGSRKPPREMLKLSGSDSSQSMNGLEVAPPGLITNFSLATAEQ
CGQETPLENMLFASFYLLDFILALVGNTLALWLFIRDHKSGTPANVFLMH
LAVADLSCVLVLPTRLVYHFSGNHWPFGEIACRLTGFLFYLNMYASIYFL
TCISADRFLAIVHPVKSLKLRRLPYAHLACAFLWVVAVAMAPLLVSPQT
VQTNHTVVCLQLYREKASHHALVSLAVAFTFPFITTVCYLLIIRSLRQG
LRVEKRLKTKAVRMIAIVLAIFLVCFVPYHVNRSVYVLHYRSHGASCATQ
RILALANRITSCLTSLNGALDPIMYFFVAEKFRHALCNLLCGKRLKGPPP
SFEGKTNESSLSAKSEL

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

GPR17

Full Name

G protein-coupled receptor 17

Alternative Names

uracil nucleotide/cysteinyl leukotriene receptor; P2Y-like receptor; R12; UDP/CysLT receptor; GPR17; G protein-coupled receptor 17

Gene ID

[2840](#)

UniProt ID

[Q13304](#)