

Product Information

MemDX™ Membrane Protein Human GPR37L1 (G protein-coupled receptor 37 like 1) Full Length

Cat. No.: **MPC0162K**

This product is for research use only and is not intended for diagnostic use.

This product is a 52.7 kDa Human GPR37L1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

GPR37L1

Protein Length

Full length

Protein Class

GPCR

Molecular Weight

52.7 kDa

TMD

7

Sequence

MRWLWPLAVSLAVILAVGLSRVSGGAPLHLGRHRAETQEQQSRSKRGTEDEAKGVQQYVPEEWAHEYPRPIHPAGLQPTKPLVATSPNPGKDGGTPDSGQELRGNLTGAPGQRLQIQNPLYPVTESSYSAYAIMLLALVVFVAVGIVGNLSVMCIVWHSYYLKSAWNLSILASLALWDFLVLFCLPIVIFNEITKQRL LGDVSCRAVPFMEVSSLGVTTFSLCALGIDRFHVATSTLPKVRPIERCQSILAKLAVIWVGSMTLAVPELLLWQLAQEPAPTMGTLDSCIMKPSASLPESLYSLVMTYQNARMWWWYFGCYFCLPILFTVTCQLVTWRVRGPPGRKSECRASHKEQCESQLNSTVVGLTVVYAFCTLPENVCNIVVAYLSTELTRQTLDLLGLINQFSTFFKGAITPVLLLCICRPLGQAFLDCCCCCCECGGASEASAANGSDNKLKTEVSSSIYFHKPRESPPLLPLGTPC

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

GPR37L1

Full Name

G protein-coupled receptor 37 like 1

Alternative Names

ETBRLP2; ETBR-LP-2; ET(B)R-LP-2; endothelin B receptor-like protein 2; endothelin type b receptor-like protein 2; prosaposin receptor GPR37L1; GPR37L1; G protein-coupled receptor 37 like 1

Gene ID

[9283](#)

UniProt ID

[O60883](#)