

Product Information

MemDX™ Membrane Protein Human GPR62 (G protein-coupled receptor 62) Full Length

Cat. No.: **MPC0172K**

This product is for research use only and is not intended for diagnostic use.

This product is a 37.6 kDa Human GPR62 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

GPR62

Protein Length

Full length

Protein Class

GPCR

Molecular Weight

37.6 kDa

TMD

7

Sequence

MANSTGLNASEVAGSLGLILAAVVEVGALLGNGALLVVLRTPGLRDALY
LAHLCVVDLLAAASIMPLGLLAAPPGLGRVRLGPAPCRAARFLSAALLP
ACTLGVAALGLARYRLIVHPLRPGSRPPVVLVTAVWAAAGLLGALSLLG
TPPAPPPAPARCSVLAGGLGPFRLWALLAFALPALLLLGAYGGIFVVAR
RAALRPPRPARGSRLHSDSLDSRLSILPPLRPRLPGGKAALAPALAVGQF
AACWLPYGCACLAPAAARAAEAAVTWVAYSFAAHPFLYGLLQRPVRLA
LGRLSRRALPGPVRCTPQAWHPRALLQCLQRPPEGPAVGPSEAPEQTPE
LAGGRSPAYQGPPSSLS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

GPR62

Full Name

G protein-coupled receptor 62

Alternative Names

GPCR8; KPG_005; G-protein coupled receptor 62; G-protein coupled receptor GPCR8; G-protein coupled receptor KPG_005; hGPCR8; probable G-protein coupled receptor 62; GPR62; G protein-coupled receptor 62

Gene ID

[118442](#)

UniProt ID

[Q9BZJ7](#)