

Product Information

MemDX™ Membrane Protein Human SLC44A2 (Solute carrier family 44 member 2) Full

Length

Cat. No.: **MPC0859K**

This product is for research use only and is not intended for diagnostic use.

This product is a 80.1 kDa Human SLC44A2 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SLC44A2

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

80.1 kDa

TMD

10

Sequence

MGDERPHYGKHGTPQKYDPTFKGPIYNRGCTDIICCVFLLAIVGYVAV
GIIAWTHGDPKVIYPTDSRGEFCGQKGTKNENKPYLFYFNIVKCASPLV
LLEFQCPTPQICVEKCPDRYLTYLNARSSRDFEYYKQFCVPGFKNNKGVA
EVLQDGDGCPAVLIPSKPLARRCFPAIHAYKGVLMVGNETTYEDGHGSRKN
ITDLVEGAKKANGVLEARQLAMRIFEDYTVSWYWIIIGLVIAMAMSLIFI
ILLRFLAGIMVWVMIMVILVLGYGIFHCYMEYSRLRGEAGSDVSLVDLG
FQTDFRVYLHLRQTLAFMIILSILEVILLIFLRKRILIAIALIKEA
SRAVGVMCSLLYPLVTFFLLCLCIAYWASTAVFLSTSNEAVYKIFDDSP
CPFTAKTCNPETFPSSNESRQCPNARCQFAFYGGESGYHRALLGLQIFNA
FMFFWLANFVLALGQVTLAGAFASYWALRKPDPLPAFPLFSAGRALRY
HTGSLAFGALILAIQIIRVILEYLDQRLKAAENKFAKCLMTCLKCCFWC
LEKFIKFLNRNAYIMIAIYGTNFCSTARNAFFLLMRNIIRVAVLDKVTDF
LFLLGKLLIVGSGILAFFFTHRIRIVQDTAPPLNYYWVPILTIVIGSY
LIAHGFFSVYGMCDVTLFLCFLEDLERNDGSAERPYPMSSTLKLLNKTN
KKA AES

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

SLC44A2

Full Name

Solute carrier family 44 member 2

Alternative Names

CTL2; PP1292; choline transporter-like protein 2; solute carrier family 44 (choline transporter), member 2; testicular tissue protein Li 47; SLC44A2; Solute carrier family 44 member 2

Gene ID

[57153](#)

UniProt ID

[Q8IWA5](#)