

Product Information

MemDX™ Membrane Protein Human SLC51B (Solute carrier family 51 subunit beta) Full Length

Cat. No.: **MPC0876K**

This product is for research use only and is not intended for diagnostic use.

This product is a 14.3 kDa Human SLC51B membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SLC51B

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

14.3 kDa

TMD

1

Sequence

MEHSEGAPGDPAGTVVPQELLEMLWFFRVEDASPWNHSILALAAVVVII
SMVLLGRSIQASRKEKMQPPEKETPEVLHLDEAKDHNSLNNLRETLLSEK
PNLAQVELELKERDVLSVFLPDVPETES

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

SLC51B

Full Name

Solute carrier family 51 subunit beta

Alternative Names

OSTB; OSTBETA; organic solute transporter subunit beta; OST-beta; organic solute transporter beta subunit; solute carrier family 51 beta subunit; SLC51B; Solute carrier family 51 subunit beta

Gene ID

[123264](#)

UniProt ID

[Q86UW2](#)