

Product Information

MemDX™ Membrane Protein Human TMEM30B (Transmembrane protein 30B) Full Length

Cat. No.: **MPC0954K**

This product is for research use only and is not intended for diagnostic use.

This product is a 38.9 kDa Human TMEM30B membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM30B

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

38.9 kDa

TMD

2

Sequence

MTWSATARGAHQPDNTAFTQQRLPAWQPILLSASIALPLFFCAGLAFIGLG
LGLYYSSNGIKELEYDYTGDPGTGNCVCAAAGQGRALPPPCSCAWYFSL
PELFQGPVYLYELTNFYQNNRRYGVSRRDQAQLSGLPSALRHPVNECAPY
QRSAAGLPIAPCGAIANSFLNDSFSLWHQRQPGGPYVEVPLDRSGIAWWT
DYHVKFRNPPLVNGSLALAFQGTAPPPNWRPVEYELSPDPNNTGFINQDF
VWWMRTAALPTFRKLYARIRQGNYSAGLPRGAYRVNITYNYPVRAFGGHK
LLIFSSISWMGGKNPFLGIAYLVVGLCILTGFVMLVVYIRYQDQDDDE
E

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM30B

Full Name

Transmembrane protein 30B

Alternative Names

CDC50B; cell cycle control protein 50B; P4-ATPase flippase complex beta subunit TMEM30B; TMEM30B; Transmembrane protein 30B

Gene ID

[161291](#)

UniProt ID

[Q3MIR4](#)