

Product Information

MemDX™ Membrane Protein Human C2CD2L (C2CD2 like) Full Length

Cat. No.: **MPC1097K**

This product is for research use only and is not intended for diagnostic use.

This product is a 76.1 kDa Human C2CD2L membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

C2CD2L

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

76.1 kDa

TMD

1

Sequence

MDPGWGQRDVGWAALLILFAASLLTVFAWLLQYARGLWLARARGDRGPGP
ALAGEPAGSLRELGVWRSLLRLRATRAGAAEPPGVRGLLASLFAFKSFRE
NWQRAWVRALNEQACRNGSSIQIAFEEVPQLPPRASISHVTCVDQSEHTM
VLRCQLSAEEVRFPVSVTQQSPAASVSMETYHVTLLPPTQLEVNLEEIPG
EGLLISWAFTDRPDLSLTVLPKLQARERGEREEQVELSTIEELIKDAIVSTQ
PAMMVNLRAC SAPGGLVPSEKPPMMPQAQPAIPRNRLFLRQLRASHLGN
ELEGTEELCCVAELDNPMQKWKTPARAGSEVEWTEDLALDLGPQSRELT
LKVLRSSSCGDTELLGQATLPVGSPSRPLSRRQLCPLTPGPGKALGPAAT
MAVELHYEEGSPRNLTPTSSTPRPSITPTKKIELDRTIMPDGTIVTTVT
TVQSRPRIDGKLDSPSRSPSKVEVTEKTTTVLSESSGPSNTSHSSSRDSH
LSNLDPVAETAIRQLTEPSGRVAKKTPTKRSTLIISGVSKVPVIAQDELA
LSLGYAASLEASVQDDAGTSGGPSSPPSDPPAMSPGPLDALSSPTSVQEA
DETTTRSDISERPSVDDIESETGSTGALETRSLKDHKVSFLRSGTKLIFRR
RPRQKEAGLSQSHDDL SNATATPSVRKKAGSFSRRLIKRFSFKSKPKANG
NPSPQL

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

C2CD2L

Full Name

C2CD2 like

Alternative Names

C2CD2L; DLNB23; TMEM24; phospholipid transfer protein C2CD2L; C2 domain-containing protein 2-like; transmembrane protein 24; KIAA0285; C2CD2 like

Gene ID

[9854](#)

UniProt ID

[O14523](#)