

Product Information

MemDX™ Membrane Protein Human PLLP (Plasmolipin) Full Length

Cat. No.: **MPC1112K**

This product is for research use only and is not intended for diagnostic use.

This product is a 19.9 kDa Human PLLP membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

PLLP

Protein Length

Full length

Protein Class

Transporter; Ion channel

Molecular Weight

19.9 kDa

TMD

4

Sequence

MAEFPSKVSTRTSSPAQGAEASVSALRPDLGFVRSRLGALMLLQLVLGLL
VWALIADTPYHLYPAYGWVMFVAVFLWLVTIVLFNLYLFQLHMKLYMVPW
PLVLMIFNISATVLYITAFIACSAAVDLTSLRGTRPYNQRAAASFFACLV
MIAYGVSAFFSYQAWRGVGSNAATSQMAGGYA

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

PLLP

Full Name

Plasmolipin

Alternative Names

PLLP; PMLP; TM4SF11; plasma membrane proteolipid (plasmolipin); transmembrane 4 superfamily member 11 (plasmolipin); Plasmolipin

Gene ID

[51090](#)

UniProt ID

[Q9Y342](#)