

Product Information

MemDX™ Membrane Protein Human TMBIM6 (Transmembrane BAX inhibitor motif containing 6) Full Length

Cat. No.: **MPC1133K**

This product is for research use only and is not intended for diagnostic use.

This product is a 26.5 kDa Human TMBIM6 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMBIM6

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

26.5 kDa

TMD

7

Sequence

MNIFDRKINF DALLKFSHITPSTQQHLKKVYASFALCMFVAAAGAYVHMV
THFIQAGLLSALGSLILMIWLMATPHSHETEQRLLAGFAFLTGVGLG
PALEFCIAVNPSILPTAFMGTA MIFTCFTLSALYARRRSYLFLGGILMSA
LSLLLLSSLGNVFFGSIWLFQANLYVGLVVMCGFVLFD TQLIIEKAEHGD
QDYIWHCIDLFLDFITVFRKLMMILAMNEKDKKKEKK

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMBIM6

Full Name

Transmembrane BAX inhibitor motif containing 6

Alternative Names

TMBIM6; BI-1; TEGT; BAXI1; bax inhibitor 1; testis enhanced gene transcript; testis-enhanced gene transcript protein; transmembrane BAX inhibitor motif-containing protein 6; Transmembrane BAX inhibitor motif containing 6

Gene ID

[7009](#)

UniProt ID

[P55061](#)