

Product Information

MemDX™ Membrane Protein Human IFITM2 (Interferon induced transmembrane protein 2)

Full Length

Cat. No.: **MPC1235K**

This product is for research use only and is not intended for diagnostic use.

This product is a 14.6 kDa Human IFITM2 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

IFITM2

Protein Length

Full length

Protein Class

Immunity

Molecular Weight

14.6 kDa

TMD

2

Sequence

MNHIVQTFSPVNSGQPPNYEMLKEEQEVAMLGVPHPNPAPPMSTVIHIRSE
TSVPDHVVWSLFTLTMNTCCLGFI AFAYSVKSRDRKMVG DVTGAQAYAS
TAKCLNIWALILGIFMTILLIIIPVLVVQAQR

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

IFITM2

Full Name

Interferon induced transmembrane protein 2

Alternative Names

IFITM2; 1-8D; DSPA2c; interferon-induced transmembrane protein 2; dispanin subfamily A member 2c; interferon-inducible protein 1-8D; Interferon induced transmembrane protein 2

Gene ID

[10581](#)

UniProt ID

[Q01629](#)