

Product Information

MemDX™ Membrane Protein Human TM9SF3 (Transmembrane 9 superfamily member 3) Full Length

Cat. No.: **MPC1244K**

This product is for research use only and is not intended for diagnostic use.

This product is a 67.8 kDa Human TM9SF3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TM9SF3

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

67.8 kDa

TMD

9

Sequence

MRPLPGALGVAAAAALWLLLLLPRTTRADEHEHTYQDKKEEVVLWMNTVGP
YHNRQETYKYFSLPFCVGSKKSISHYHETLGEALQGVELEFSGLDIKFKD
DVMPATYCEIDLDKEKRDAFVYAIAKNHYWYQMYIDDLPIWGIVGEADENG
EDYYLWTYKKLEIGFNGNRIVDVNLTSEGKVKLPNTKIQMSYSVKWKKS
DVKFEDRFDKYLDPSFFQHRIHWFSIFNSFMMVIFLVGLVSMILMRTLK
DYARYSKEEEMDDMDRDLGDEYGWKQVHGDVFRPSSHPLIFSSLIGSGCQ
IFAVSLVIVAMIEDLYTERGSMLSTAIFVYAATSPVNGYFGGSLYARQ
GGRRWIKQMFIGAFLIPAMVCGTAFFINFIAIYYHASRAIPFGTMVAVCC
ICFFVILPLNLVGTILGRNLGQPNFPCRVNAVPRPIPEKKWFMEPAVIV
CLGGILPFGSIFIEMFYIFTSFWAYKIYYVYGFMMLVLVILCIVTVCVTI
VCTYFLLNAEDYRWQWTSFLSAASTAIYVYMYSFYFFKTKMYGLFQTS
FYFGYMAVFSTALGIMCGAIGYMGTSAFVRKIYTNVKID

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TM9SF3

Full Name

Transmembrane 9 superfamily member 3

Alternative Names

TM9SF3; SMBP; EP70-P-iso; SM-11044-binding protein; dinucleotide oxidase disulfide thiol exchanger 3 superfamily member 3; endomembrane protein emp70 precursor isolog; Transmembrane 9 superfamily member 3

Gene ID

[56889](#)

UniProt ID

[Q9HD45](#)