

Product Information

MemDX™ Membrane Protein Human TMED2 (Transmembrane p24 trafficking protein 2) Full Length

Cat. No.: **MPC1259K**

This product is for research use only and is not intended for diagnostic use.

This product is a 22.7 kDa Human TMED2 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMED2

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

22.7 kDa

TMD

1

Sequence

MVTLAELLVLLAALLATVSGYFVSIDAHAECCFFERVTSGTKMGLIFEVA
EGGFLDIDVEITGPDNKGIIYKGDRESSGKYTFAAHMDGTYKFCFSNRMST
MTPKIVMFTIDIGEAPKGQDMETEAHQNKLEEMINELAVAMTAVKHEQEY
MEVRERIHRAINNTNSRVVLWSFFEALVLVAMTLGQIYYLKRFFEVRV
V

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMED2

Full Name

Transmembrane p24 trafficking protein 2

Alternative Names

TMED2; p24; P24A; RNP24; p24b1; p24beta1; transmembrane emp24 domain-containing protein 2; coated vesicle membrane protein; membrane protein p24A; p24 family protein beta-1; transmembrane emp24 domain trafficking protein 2; Transmembrane p24 trafficking protein 2

Gene ID

[10959](#)

UniProt ID

[Q15363](#)