

# Product Information

## MemDX™ Membrane Protein Human TMEM201 (Transmembrane protein 201) Full Length

Cat. No.: **MPC1260K**

This product is for research use only and is not intended for diagnostic use.

This product is a 72.2 kDa Human TMEM201 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

TMEM201

#### Protein Length

Full length

#### Protein Class

Transporter

#### Molecular Weight

72.2 kDa

#### TMD

5

#### Sequence

MEGVSALLARCPTAGLAGGLGVTACAAAGVLLYRIARRMKPTHTMVNCWF  
CNQDTLVPYGNRNCWDCPHCEQYNGFQENG DY NKPIPAQYLEHLNHVVSS  
APSLRDPSQPQQWVSSQVLLCKRCNHHQTTKIKQLAAAFAPREEGRYDEEV  
EVYRHHLEQMYKLCRPCQAAVEYYIKHQNRQLRALLSHQFKRREADQTH  
AQNFSSAVKSPVQVILLRALAFLACAFLLTTALYGASGHFAPGTTVPLAL  
PPGGNGSATPDNGTTPGAEGWRQLLGLLPEHMAEKLCEAWAFGQSHQTGV  
VALGLLTCLLAMLLAGRI RLRRIDAFCTCLWALLGLHLAEQHLQAASPS  
WLDTLKFSTTSLCCLVGFTA AVATRKATGPRRFRPRRFFPGDSAGLFPTS  
PSLAIPHPVGGSPASLFIPSPSFLPLANQQLFRSPRRTSPSSLPGRLS  
RALSLGTIPSLTRADSGYLFSGSRPPSQVSRSGEFPVSDYFSLLSGSCPS  
SPLPSPAPSVAGSVASSSGSLRHRRPLISPARLNLKGQKLLFPSPGGEA  
PTTPSSSDEHSPHNGSLFTMEPPHVPRKPPLQDVKHALDLRSKLERGSAC  
SNRSIKKEDDSSQSSTCVVDTTTRGCSEEAATWRGRFGPSLVRGLLAVSL  
AANALFTSVFLYQSLR

### Product Description

#### Expression Systems

HEK293

**Tag**

Based on specific requirements

**Protein Format**

Detergent or based on specific requirements

**Form**

Liquid

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

TMEM201

**Full Name**

Transmembrane protein 201

**Alternative Names**

TMEM201; Ima1; NET5; SAMP1; RP13-15M17.2; spindle-associated membrane protein 1; Transmembrane protein 201

**Gene ID**

[199953](#)

**UniProt ID**

[Q5SNT2](#)