

Product Information

MemDX™ Membrane Protein Human TMED5 (Transmembrane p24 trafficking protein 5) Full Length

Cat. No.: **MPC1276K**

This product is for research use only and is not intended for diagnostic use.

This product is a 26 kDa Human TMED5 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMED5

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

26 kDa

TMD

1

Sequence

MGDKIWLPFPVLLLAALPPVLLPGAAGFTPSLDSDFTFPLPAGQKECFYQ
PMPLKASLEIEYQVLDGAGLDIDFHLASPEGKTLVFEQRKSDGVHTVETE
VGDMFCFDNTFSTISEKVIFFELIDNMGEQAQEQEDWKKYITGTDILD
MKLEDILESINSIKSRLSKSGHIQTLLRAFEARDRNIQESNFD RVNFWSM
VNLVVMVVVSAIQVYMLKSLFEDKRKSRT

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMED5

Full Name

Transmembrane p24 trafficking protein 5

Alternative Names

TMED5; p28; p24g2; CGI-100; transmembrane emp24 domain-containing protein 5; p24 family protein gamma-2; p24gamma2; transmembrane emp24 protein transport domain containing 5; Transmembrane p24 trafficking protein 5

Gene ID

[50999](#)

UniProt ID

[Q9Y3A6](#)