

Product Information

MemDX™ Membrane Protein Human TMED3 (Transmembrane p24 trafficking protein 3) Full Length

Cat. No.: **MPC1282K**

This product is for research use only and is not intended for diagnostic use.

This product is a 24.7 kDa Human TMED3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMED3

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

24.7 kDa

TMD

1

Sequence

MGSTVPRASVLLLLLLRRAEQPCGAELTFELPDNAKQCFHEEVEQGKVK
FSLDYQVITGGHYDVCYVEDPQGNTIYRETKKQYDSFTYRAEVKGVYQF
CFSNEFSTFSHKTVYFDFQVGDEPPILPDMGNRVLTALTQMESACVTIHEA
LKTVIDSQTHYRLREAQDRARAEDLNSRVSYWSVGETIALFVVSFSQVLL
LKSFFTEKRPISRAVHS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMED3

Full Name

Transmembrane p24 trafficking protein 3

Alternative Names

TMED3; p26; P24B; p24g4; C15orf22; p24gamma4; transmembrane emp24 domain-containing protein 3; integral type I protein; membrane protein p24B; p24 family protein gamma-4; testis tissue sperm-binding protein Li 48e; transmembrane emp24 protein transport domain containing 3; Transmembrane p24 trafficking protein 3

Gene ID

[23423](#)

UniProt ID

[Q9Y3Q3](#)