

Product Information

MemDX™ Membrane Protein Human TMEM117 (Transmembrane protein 117) Full Length

Cat. No.: **MPC1351K**

This product is for research use only and is not intended for diagnostic use.

This product is a 60.1 kDa Human TMEM117 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM117

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

60.1 kDa

TMD

8

Sequence

MGKDFRYYFQHPWSRMIVAYLVIFFNFLIFAEDPVSHSQTEANVIVVGNC
FSFVTNKYPRGVGWRILKVLLWLLAILTGLIAGKFLFHQRLFGQLRLKLM
FREDHGSWMTMFFSTILFLFIFSHIYNTILLMDGNMGAYIITDYMGIIRNE
SFMKLAAGVTWMGDFVTAWMVTDMMMLQDKPYPDWGKSARAFWKKGNVRIT
LFWTVLFTLTSSVVVLVITTDWISWDKLNRGFLPSDEVSRFLASFILVFD
LLIVMQDWEFPHFMGDVDVNLPGLHTPHMQFKIPFFQKIFKEEYRIHITG
KWFNYGIIFLVLILDNLNMWKNQIFYKPHEYGQYIGPGQKIYTVKDSESLK
DLNRTKLSWEWRSNHTNPRTNKTYVEGDMFLHSRFIGASLDVKCLAFVPS
LIAFVWFGFFIWWFFGRFLKNEPRMENQDKTYTRMKRKSPSEHSKDMGTR
ENTQASVEDPLNDPSLVCIRSDFNEIVYKSSHLTSENLSQLNESTSATE
ADQDPTTSKSTPTN

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM117

Full Name

Transmembrane protein 117

Alternative Names

TMEM117; Transmembrane protein 117

Gene ID

[84216](#)

UniProt ID

[Q9H0C3](#)