

Product Information

MemDX™ Membrane Protein Human TMBIM4 (Transmembrane BAX inhibitor motif containing 4) Full Length

Cat. No.: **MPC1356K**

This product is for research use only and is not intended for diagnostic use.

This product is a 26.9 kDa Human TMBIM4 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMBIM4

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

26.9 kDa

TMD

7

Sequence

MADPDPRYPRSSIEDDFNYGSSVASATVHIRMAFLRKVYSILSLQVLLTT
VTSTVFLYFESVRTFVHESPALILLFALGSLGLIFALILNRHKYPLNLYL
LFGFTLLEALTAVVVTFYDVYIILQAFILTTTVFFGLTVYTLQSKKDFS
KFGAGLFALLWILCLSGFLKFFFYSEIMELVAAAAGALLFCGFIYDTHS
LMHKLSPEEYVLA AISLYLDIINLFLHLLRFLEAVNKK

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMBIM4

Full Name

Transmembrane BAX inhibitor motif containing 4

Alternative Names

TMBIM4; S1R; GAAP; LFG4; ZPRO; CGI-119; protein lifeguard 4; Golgi anti-apoptotic protein; Z-protein; transmembrane BAX inhibitor motif-containing protein 4; Transmembrane BAX inhibitor motif containing 4

Gene ID

[51643](#)

UniProt ID

[Q9HC24](#)