

Product Information

MemDX™ Membrane Protein Human TMEM184C (Transmembrane protein 184C) Full Length

Cat. No.: **MPC1380K**

This product is for research use only and is not intended for diagnostic use.

This product is a 50.1 kDa Human TMEM184C membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM184C

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

50.1 kDa

TMD

7

Sequence

MPCTCTWRNWRQWIRPLVAVIYLVSIIVAVPLCVWELQKLEVGIHTKAWF
IAGIFLLLTIPISLWVILQHLVHYTQPELQKPIIRILWMVPIYSLDSWIA
LKYPGIAIYVDTCRECYEAYVIYNFMGFLTNYLTNRYPNLVLILEAKDQQ
KHFPPLCCPPWAMGEVLLFRCKLGVLYTVVRPFTTIVALICELGIIYD
EGNFSFSNAWTYLVIINNMSQLFAMYCLLLFYKVLKEELSPIQPVGKFLC
VKLVVFSFWQAVVIALLVKVGVISEKHTWEWQTVEAVATGLQDFIICIE
MFLAAIAHHYTF SYKPYVQEAEEGSCFDSFLAMWDVSDIRDDISEQVRHV
GRTVRGHPRKLPEDQDQNEHTSLLSSSSQDAISIASSMPPSPMGHYQG
FGHTVTPQTTPPTAKISDEILSDTIGEKKEPSDKSVDS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM184C

Full Name

Transmembrane protein 184C

Alternative Names

TMEM184C; TMEM34; transmembrane protein 34; Transmembrane protein 184C

Gene ID

[55751](#)

UniProt ID

[Q9NVA4](#)