

# Product Information

## MemDX™ Membrane Protein Human TMEM192 (Transmembrane protein 192) Full Length

Cat. No.: **MPC1384K**

This product is for research use only and is not intended for diagnostic use.

This product is a 30.9 kDa Human TMEM192 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

TMEM192

#### Protein Length

Full length

#### Protein Class

Transporter

#### Molecular Weight

30.9 kDa

#### TMD

4

#### Sequence

MAAGGRMEDGSLDITQSIEDDPLLDAQLLPHHSLQAHFRPRFHPLPTVII  
VNLLWFIHLVVFVLAFLTGVLCSYPNPNEKCPGNYTNPLKVQTVIILGK  
VILWILHLLLECYIQYHHSKIRNRGYNLIYRSTRHLKRLALMIQSSGNTV  
LLLILCMQHSFPEPGRLYLDLILAILALELICSLICLLIYTVKIRRFNKA  
KPEPDILEEEKIYAYPSNITSETGFRTISSLEEIVEKQGDITIEYLKRHNA  
LLSKRLLALTSSDLGCQPSRT

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

#### Protein Format

Detergent or based on specific requirements

**Form**

Liquid

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

TMEM192

**Full Name**

Transmembrane protein 192

**Alternative Names**

TMEM192; Transmembrane protein 192

**Gene ID**

[201931](#)

**UniProt ID**

[Q8IY95](#)