

Product Information

MemDX™ Membrane Protein Human LRIT2 (Leucine rich repeat, Ig-like and transmembrane domains 2) Full Length

Cat. No.: **MPC1420K**

This product is for research use only and is not intended for diagnostic use.

This product is a 60.1 kDa Human LRIT2 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

LRIT2

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

60.1 kDa

TMD

1

Sequence

MASVFHYFLLVLVFLDTHAAQPFCLPGCTCSEESFGRTLQCTSVSLGKIP
GNLSEEFKQVRIENSPLFEMPQGSFINMSTLEYLWLNFNINISVIHLGALE
HLPRLRELRLLEGKLCSPWTAFRATPLLRVLDLKRNKIDALPELALQFL
VSLTYLDLSSNRLTVVSKSVFLNWPAYQKCRQPDCAEILSSLVVALHDN
PWVCDCLRLGLVQFVKISITLPVILVNSYLICQGPLSKAGQLFHETELSAC
MKPQISTPSANITIRAGQNVTLRCLAQASPSPSIAWTYPLSMWREFDVL
TSTGEDTALSELAIPTAAHLVDSGNYTCMASNSIGKSNLVISLHVQPAQAL
HAPDLSIPSEGNAYIDLRVVKQTVHGILLEWLAVADTSKEEWFTLYIAS
DEAFRKEVVHIGPGINTYAVDDLLPGTKYEACLSLEGQPPHQQGQCVAFVT
GRDAGGLEAREHLLHVTVLLCVLLAVPVGAYAWAAQGPCSCSKWVLRGC
LHRRKAPSCTPAAPQSKDGSFREHPAVCDDGEGHIDTEGDKEKGGTEDNS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

LRIT2

Full Name

Leucine rich repeat, Ig-like and transmembrane domains 2

Introduction

LRRC22; leucine-rich repeat, immunoglobulin-like domain and transmembrane domain-containing protein 2; leucine rich repeat containing 22; leucine-rich repeat, immunoglobulin-like and transmembrane domains 2; leucine-rich repeat-containing protein 22

Alternative Names

LRIT2; Leucine rich repeat, Ig-like and transmembrane domains 2

Gene ID

[340745](#)

UniProt ID

[A6NDA9](#)