

Product Information

MemDX™ Membrane Protein Human TMEM11 (Transmembrane protein 11) Full Length

Cat. No.: **MPC1424K**

This product is for research use only and is not intended for diagnostic use.

This product is a 21.5 kDa Human TMEM11 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM11

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

21.5 kDa

TMD

2

Sequence

MAAWGRRRLGPGSSGGSARERVSLSATDCYIVHEIYNGENAQDQFEYELE
QALEAQYKYIVIEPTRIGDETARWITVGNCLHKTAVLAGTACLFPLALP
LDYSHYISLPAGVLSLACCTLYGISWQFDPCKYQVEYDAYKLSRLPLHT
LTSSTPVVLVRKDDLHRKRLHNTIALAALVYCVKKIYELYAV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM11

Full Name

Transmembrane protein 11

Alternative Names

TMEM11; PM1; PMI; C17orf35; transmembrane protein 11, mitochondrial; putative receptor protein; Transmembrane protein 11

Gene ID

[8834](#)

UniProt ID

[P17152](#)