

Product Information

MemDX™ Membrane Protein Human CD19 (CD19 molecule) Full Length

Cat. No.: **MPC1432K**

This product is for research use only and is not intended for diagnostic use.

This product is a 61.1 kDa Human CD19 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

CD19

Protein Length

Full length

Protein Class

Immunity

Molecular Weight

61.1 kDa

TMD

1

Sequence

MPPPRLLFLLFLTPMEVRPEEPLVVKVEEGDNAVLQCLKGTSDGPTQQL
TWSRESPLKPFLLKSLGLPGLGIHMRPLAIWLFIFNVSQQMGGFYLCQPG
PPSEKAWQPGWTVNVEGSGELFRWNVSDLGGLGCLKNRSSEGPSSPSGK
LMSPKLYVWAKDRPEIWEGEPPCLPPRDSLNLQSLSQDLTMAPGSTLWLS
GVPPDSVSRGPLSWTHVHPKGPKSLLSLELKDDRPARDMWVWVETGLLLPR
ATAQDAGKYYCHRGNLMSFHLEITARPVLWHWLLRTGGWKVSAVTLAYL
IFCLCSLVGILHLQRALVLRKRKRMTDPTRRFFKVTPPPGSGPQNQYGN
VLSLPTPTSGLGRAQRWAAGLGGTAPSYGNPSSDVQADGALGSRSPPGVG
PEEEEGEGYEEPDSSEDESEFYENDSNLGQDQLSQDGSYENPEDEPLGPE
DEDSFSNAESYENEDEELTQPVARTMDFLSPHGSADWPSREATSLGSQSY
EDMRGILYAAPQLRSIRGQPGPNHEEDADSYENMDNPDGPDPAWGGGGRM
GTWSTR

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

CD19

Full Name

CD19 molecule

Introduction

This gene encodes a member of the immunoglobulin gene superfamily. Expression of this cell surface protein is restricted to B cell lymphocytes. This protein is a reliable marker for pre-B cells but its expression diminishes during terminal B cell differentiation in antibody secreting plasma cells. The protein has two N-terminal extracellular Ig-like domains separated by a non-Ig-like domain, a hydrophobic transmembrane domain, and a large C-terminal cytoplasmic domain. This protein forms a complex with several membrane proteins including complement receptor type 2 (CD21) and tetraspanin (CD81) and this complex reduces the threshold for antigen-initiated B cell activation. Activation of this B-cell antigen receptor complex activates the phosphatidylinositol 3-kinase signalling pathway and the subsequent release of intracellular stores of calcium ions. This protein is a target of chimeric antigen receptor (CAR) T-cells used in the treatment of lymphoblastic leukemia. Mutations in this gene are associated with the disease common variable immunodeficiency 3 (CVID3) which results in a failure of B-cell differentiation and impaired secretion of immunoglobulins. CVID3 is characterized by hypogammaglobulinemia, an inability to mount an antibody response to antigen, and recurrent bacterial infections. Alternative splicing results in multiple transcript variants encoding distinct isoforms.

Alternative Names

CD19; B4; CVID3; B-lymphocyte antigen CD19; B-lymphocyte surface antigen B4; T-cell surface antigen Leu-12; differentiation antigen CD19; CD19 molecule

Gene ID

[930](#)

UniProt ID

[P15391](#)